

# PCR Diagnostic Tests for Companion Animals

## CANINE

- Babesia Gibsoni* RM 90.10
- \**Leptospira* spp. RM 90.10
- \*Canine Parvovirus RM 90.10

## FELINE

- Feline Immunodeficiency Virus (FIV) RM 106.00
- \*Feline Leukemia Virus (FeLV) RM 106.00
- \*Feline Parvovirus RM 90.10

Pets are part of our family, we love them and we try our best to protect them from any source of illness. In line with that, VFAD presents our new PCR (Polymerase Chain Reaction) diagnostic tests for companion animals to provide accurate detection and enable pet owners to seek early treatment.

## FEATURES




Accurate Result



Cost-Effective



Fast (within 24 hours)



**TECHNOLOGY**  
**Polymerase Chain Reaction (PCR)**

## SAMPLE SUBMISSION

- Babesia gibsoni* : 1 ml of whole blood
- Leptospira* spp. : 1 ml of whole blood (first 7-10 days of infection), 10 ml of urine.
- Canine Parvovirus : Minimum 2 grams of fresh faecal samples
- Feline Immunodeficiency Virus : 1 ml of whole blood
- Feline Leukemia Virus : 1 ml of whole blood
- Feline Parvovirus : Minimum 2 grams of fresh faecal samples

Note: Blood and faecal samples should be stored and transported at 4 °C.  
(\* ) is not included in the SAMM Accreditation Schedule of this laboratory.



Vet Food Agro Diagnostics (M) Sdn Bhd  
200101024963 (560721-X)

Lot 18B, Jalan 241, Seksyen 51A, 46100 Petaling Jaya, Selangor, Malaysia  
(T) 603 7873 6405/7873 7355 (372/379) (F) 603 2770 0119  
(W) www.vfad.com.my (E) vfad@rhonema.com



Get your lab results online at VETLAB® Reporting System (VRS) Online: [www.vfad.com.my/vrsonline](http://www.vfad.com.my/vrsonline)